


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Methylphenol - analytical standard

Product Number: 442391
CAS Number: 108-39-4
Formula: C₇H₈O
Formula Weight: 108.14 g/mol



TEST

Specification

Infrared spectrum
Infrared spectrum
Purity (GC)
GC - Mass Spec
Expiration Date Period
3 Years

Conforms to Structure
Conforms to Structure
≥ 96.0 %
Matches: NIST

Specification: PRD.0.ZQ5.10000048035

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.