

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
KT 5823 - ≥85% (HPLC), lyophilized powder

Product Number: K1388
CAS Number: 126643-37-6
MDL: MFCD09878278
Formula: C₂₉H₂₅N₃O₅
Formula Weight: 495.53 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Solid
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
5 mg/mL, Ethyl Acetate	
Purity (HPLC)	≥ 85 %

Specification: PRD.1.ZQ5.10000032266

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.