

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methyl cis-10-heptadecenoate - ≥99% (GC), liquid

Product Number: H9021
CAS Number: 75190-82-8
MDL: MFCD00133176
Formula: C₁₈H₃₄O₂
Formula Weight: 282.46 g/mol
Storage Temperature: -20 °C



TEST	Specification
Proton NMR spectrum	Conforms to Structure
Purity (GC)	≥ 99 %
Recommended Retest Period 2 years	-----
Appearance (Color)	Colorless
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid

Specification: PRD.0.ZQ5.10000059415

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.