

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

(-)-Isopinocampheylborane TMEDA complex

Product Number:

251496

CAS Number:

67826-92-0

MDL:


MFCD00077801

Formula:

C₁₃H₂₇BN

Formula Weight:

208.17 g/mol



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Crystals to Crystalline Powder to Powder

Infrared spectrum

Conforms to Structure

Carbon

73.5- 75.5 %

Nitrogen

6.0- 7.0 %

Optical Rotation

65.0-- 55.0-

Degrees(C= 1%, CHCL3)

Remarks:

Specification Date : 06-05-2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.