

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

3-Chloro-3-methyl-1-butyne 97%

Product Number:

301345

CAS Number:

1111-97-3

MDL:

MFCD00190221

Formula:


C₅H₇Cl

Formula Weight:

102.56 g/mol

Storage Temperature:

2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	96.5- 103.5 %
Purity (GC)	≥96.5 %

Remarks:

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.