


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Amino-1,3-dimethylpyrazole - 97%

Product Number: 532223
CAS Number: 3524-32-1
MDL: MFCD00051651
Formula: C₅H₉N₃
Formula Weight: 111.15 g/mol



TEST

Specification

Appearance (Color)
Yellow to Orange to Tan
Appearance (Form)
Infrared spectrum
Purity (Titration by HCL04)
Purity (GC)

Conforms to Requirements
Powder or Crystals
Conforms to Structure
96.5 - 103.5 %
≥ 96.5 %

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.