

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

1.3-Dimethyl-8-phenylxanthine crystalline

Product Number:

P2278 961-45-5

CAS Number: MDL:

MFCD00005582 C13H12N4O2

Formula: Formula Weight:

256.26 g/mol

H.C. N N N N N N N N N N N N N N N N N N	
	$\overline{}$

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	59.4- 62.4 %
Nitrogen	21.3- 22.4 %
Purity (HPLC)	≥ 98.0 %
Solubility (Turbidity) c = 50 mg/ml; 1N NaOH:Ethanol (1:3)	Clear

Remarks:

Solubility (Color)

Specification Date: 10/04/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Yellow