


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,5-Dihydroxybenzaldehyde - 98%

Product Number: D108200
CAS Number: 1194-98-5
MDL: MFCD00003333
Formula: C₇H₆O₃
Formula Weight: 138.12 g/mol



TEST

Specification

Appearance (Color)

Faint Yellow to Dark Yellow and
Yellow-Green and Green-Yellow

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Titration by Hydroxylamine

97.5 - 102.5 %

Specification: PRD.1.ZQ5.10000026424

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.