

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Maltohexaose - $\geq 65\%$ (HPLC)

Product Number: M9153
CAS Number: 34620-77-4
MDL: MFCD00083650
Formula: C₃₆H₆₂O₃₁
Formula Weight: 990.86 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	$\leq 8\%$
Purity (HPLC)	$\geq 65\%$

Specification: PRD.2.ZQ5.10000008036

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.