

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,5-Dimethylphenylboronic acid

Product Number: 483516
CAS Number: 85199-06-0
MDL: MFCD01863525
Formula: C₈H₁₁BO₂
Formula Weight: 149.98 g/mol

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Proton NMR Spectrum	Conforms to Structure
Purity (Titration by NaOH)	≥ 95.0 %

Specification: PRD.0.ZQ5.10000039381

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.