


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isophthalic acid - 99%

Product Number: I19209
CAS Number: 121-91-5
MDL: MFCD00002516
Formula: C₈H₆O₄
Formula Weight: 166.13 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Titration with NaOH	98.5 - 101.5 %
Purity (GC)	≥ 98.5 %
Solubility (Turbidity) 50mg/ml, 1N NAOH	Clear
Solubility (Color)	Colorless to Faint Yellow

Specification: PRD.1.ZQ5.10000008511

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.