

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,3-Cyclopentanediol, mixture of cis and trans - 95%

Product Number:

192805

CAS Number:

59719-74-3

MDL:


MFCD00001367

Formula:

C₅H₁₀O₂

Formula Weight:

102.13 g/mol



TEST

Specification

Appearance (Color)

Colorless

Appearance (Form)

Confirmed

Viscous Liquid, Semi-Solid, or Solid

Infrared spectrum

Conforms to Structure

Refractive index at 20 °C

1.481 - 1.487

GC Purity (sum of isomers)

≥ 94.0 %A

GC (% Cis Isomer)

GC (% Trans Isomer)

Specification: PRD.0.ZQ5.10000036069

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.