

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sulfamic acid - ReagentPlus®, ≥99%

Product Number: 242772
CAS Number: 5329-14-6
MDL: MFCD00011603

Formula Weight: 97,09 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Titration with NaOH	99.0 - 101.0 %
Solubility (Turbidity) 50mg/ml, H2O	Clear
Solubility (Color)	Colorless
Recommended Retest Period 6 Years	-----

Specification: PRD.0.ZQ5.10000095387

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.