

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(tetrafluoroethylene) – powder, 35 µm particle size

Product Number:

468096


CAS Number:

9002-84-0

MDL:

MFCD00084455

Formula:

(CF₂CF₂)_n

TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Mean Particle Diameter

31.00 - 36.00 micron

Particle Size

≤ 160.00

99% of Total, micron

Specification: PRD.0.ZQ5.10000026440

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.