

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Okadaic acid sodium salt from *Prorocentrum concavum* - ≥90% (HPLC), film

Product Number:

O7760

CAS Number:

209266-80-8

MDL:

MFCD00210208

Formula:


C₄₄H₆₇NaO₁₃

Formula Weight:

826.98 g/mol

Storage Temperature:

-20 °C



TEST	Specification
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
250 ug/mL, 0.1% HOAc in CH ₃ CN:H ₂ O (1:1)	
Identity by ICP - Sodium	Positive
Purity (HPLC)	≥ 90 %

Specification: PRD.1.ZQ5.10000036930

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.