


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4,4'-Diphenyl-2,2'-dipyridyl - technical grade

Product Number: 366323
CAS Number: 6153-92-0
MDL: MFCD00006421
Formula: C₂₂H₁₆N₂
Formula Weight: 308.38 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Carbon
Nitrogen

Yellow to Tan
Powder
Conforms to Structure

Specification: PRD.0.ZQ5.1000020809

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.