

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium sulfate - ≥99.0%, suitable for plant cell culture

Product Number: S5640
CAS Number: 7757-82-6
MDL: MFCD00003504
Formula: Na₂O₄S
Formula Weight: 142.04 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
10 g plus 100 ml H ₂ O	
Loss on Drying	≤ 0.5 %
% Purity (Titration)	≥ 99.0
Plant Cell Culture Test	Pass

Specification: PRD.1.ZQ5.10000036798

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

