


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium benzoate - ReagentPlus®, 99%

Product Number: 109169
CAS Number: 532-32-1
MDL: MFCD00012463
Formula: C₇H₅NaO₂
Formula Weight: 144.10 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder or Crystals or Granules
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	98.5 - 101.5 %
Purity (HPLC)	≥ 98.5 %
Recommended Retest Period	-----
4 Years	

Specification: PRD.1.ZQ5.10000003844

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

