

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N,N'-Bis(benzyloxycarbonyl)-1H-pyrazole-1-carboxamide - 97%

Product Number:

439428

CAS Number:

152120-55-3

MDL:


MFCD02683516

Formula:

C₂₀H₁₈N₄O₄

Formula Weight:

378.38 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals, Crystalline Powder and/or Chunks

Infrared spectrum

Conforms to Structure

Carbon

61.3 - 65.7 %

Nitrogen

14.3 - 15.3 %

Purity (HPLC)

≥ 96.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.