

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N-Boc-1H-pyrazole-4-boronic acid

Product Number: **715263**
CAS Number: 947533-31-5

Formula: C₈H₁₃BN₂O₄
Formula Weight: 212.01 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color) White to Beige to Yellow	Conforms to Requirements
Appearance (Form) Powder, Crystals or Crystalline Powder and/or Chunks	Conforms to Requirements
Proton NMR Spectrum May Contain Varying Amounts of Anhydride	Conforms to Structure
Purity (Titration by NaOH)	≥ 95.0 %

Specification: PRD.0.ZQ5.10000039387

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.