

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-Human IgG (Fab specific)–FITC antibody produced in goat – affinity isolated antibody, buffered aqueous solution
Product Number: F5512

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	3.0 - 6.5 mg/ml
F/P Molar Ratio	2.5 - 6.5
IF Cells	Hep-2 Human Cells
IF-Titer	1:16-1:32
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000015297

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.