


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2-Naphthaleneethanol - 98%

Product Number: **188107**  
CAS Number: 1485-07-0  
MDL: MFCD00004128  
Formula: C<sub>12</sub>H<sub>12</sub>O  
Formula Weight: 172.22 g/mol



### TEST

### Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Initial Melting Point	
Specification: Report Only	
Final Melting Point	
Specification: Report Only	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %

Specification Date : 04/27/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.