

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Echistatin from Echi carinatus – lyophilized powder, γ -irradiated, BioXtra

Product Number:

E1518

CAS Number:

129038-42-2

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White
Solubility (Solvent)	Water
Solubility (Conc)	0.10 - 0.11 mg/ml
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
Gel Electroph. (Gel)	SDS-Page
Gel Electroph. (Stain)	Coomassie
Gel Electroph. (MW Major Band)	5.0 - 6.0 kDa
Gel Electroph. (Purity)	\geq 95.00 %
Gel Electroph. (Description)	1 Major Band
Vial Content (Size)	50.00 - 50.00
Vial Content (Actual Content)	\geq 50.00
Vial Content (Unit Of Measure)	μ g Protein
Vial Content (Testing Method)	Lowry
Gamma-Irrad.(Sterilization)	Pass
Gamma-Irrad.(Measured Dose)	Min. 2.5 Mrad (25KGy)
Recommended Retest Period	-----
3 Years	

Specification: PRE.2.ZQ5.10000000832

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.