

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Lactoferricin B, fragment 4-14 trifluoroacetate salt - ≥95% (HPLC)

Product Number:

L1290

MDL:

MFCD04040902

Formula:

C70H113N25O13S

Formula Weight:

1,544.87 g/mol

Storage Temperature:

-20 °C

| TEST | Specification |
|------|---------------|
|------|---------------|

Appearance (Color)

White

Appearance (Form)

Powder

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

1 mg/ml, H<sub>2</sub>O

% Purity

≥ 95

% Peptide content

≥ 55

Specification: PRD.0.ZQ5.10000037428

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.