

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

trans-2-Hexen-1-al - natural, ≥95%, FG

Product Number:

W256110

CAS Number:

6728-26-3

MDL:

MFCD00007008

Formula:


C₆H₁₀O

Formula Weight:

98.14 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Refractive index at 20 °C

1.445 - 1.449

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 95.00 %

Arsenic (As)

≤ 3 ppm

Cadmium (Cd)

≤ 1 ppm

Mercury (Hg)

≤ 1 ppm

Lead (Pb)

≤ 10 ppm

Expiration Date Period

5 Years

Specification: PRD.0.ZQ5.10000044676

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.