

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium nitrite - 99.999% trace metals basis

Product Number: 563218
CAS Number: 7632-00-0
MDL: MFCD00011118
Formula: NNaO₂
Formula Weight: 69.00 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals
Titration by KMNO ₄	≥ 98.0 %
ICP Major Analysis	Confirmed
Confirms Sodium Component	
Purity	Conforms
99.999% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.1.ZQ5.10000004855

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.