

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Butyl vinyl ether – contains 0.01% potassium hydroxide as stabilizer, 98%

Product Number: 110299
CAS Number: 111-34-2
MDL: MFCD00009454
Formula: C₆H₁₂O
Formula Weight: 100.16 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %
Stabilizer ~ 100ppm KOH	Confirmed

Specification: PRD.0.ZQ5.10000042223

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.