


Product Specification

Product Name:
3,3'-Diaminobenzidine - 98%

Product Number: D12384
CAS Number: 91-95-2
Formula: C₁₂H₁₄N₄
Formula Weight: 214,27 g/mol



TEST	Specification
------	---------------

Appearance (Colour)	Off White to Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Purity (HPLC)	≥ 97.5 %

Specification: PRE.3.ZQ5.10000021252

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

