

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Copper(II) chloride - anhydrous, powder, ≥99.995% trace metals basis

Product Number: 451665
CAS Number: 7447-39-4
MDL: MFCD00010972
Formula: Cl₂Cu
Formula Weight: 134.45 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Brown
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Cu Component	
Trace Metal Analysis	≤ 50.0 ppm
Purity	Conforms
> = 99.995% Based On Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000023450

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.