

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Starch, soluble - ACS reagent

Product Number: **S9765**
 CAS Number: 9005-84-9
 MDL: MFCD00082026
 Formula: (C₆H₁₀O₅)_n
 Formula Weight: 342.30 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility	Pass
Note	-----
ACS tests	
pH	5.0 - 7.0
2% solution at 25 deg C	
Residue on ignition (Ash)	≤ 0.4 %
Sensitivity Test	Pass
Recommended Retest Period	-----
2 years	

Specification: PRD.2.ZQ5.10000010928

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.