


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Chloroanthraquinone - 97%

Product Number: 156299
CAS Number: 131-09-9
MDL: MFCD00001227
Formula: C₁₄H₇ClO₂
Formula Weight: 242.66 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Grey-Green to Grey-Yellow, Yellow-Tan or Grey-Brown

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

96.5 - 103.5 %

After Oxygen Combustion

Specification: PRD.0.ZQ5.10000042265

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.