

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Nitric Oxide Synthase, Universal antibody produced in mouse - clone NOS-3F7-B11-B5, ascites fluid

Product Number: N218

MDL: MFCD01324876

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Ascites fluid containing < or = 0.1% Sodium Azide

Working Dilution	Conforms
------------------	----------

by Western Blot

1:50 - 1:500

Working Dilution	Conforms
------------------	----------

by Immunohistochemistry

1:100

Specification: PRD.1.ZQ5.10000046335

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.