

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
GAPDH - standard for protein electrophoresis

Product Number: **G5262**  
CAS Number: 9001-50-7  
MDL: MFCD00131203

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
mg protein/vial (Biuret)	5 - 7
Capillary Electrophoresis	Suitable
Suitable as an electrophoresis marker. One well defined major band with MW of 36000(+ /- 15%).	

Specification: PRD.0.ZQ5.10000076432

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.