

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Bovine IgG (whole molecule)-FITC antibody produced in rabbit - affinity isolated antibody, buffered aqueous solution

Product Number:

F7887

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	3.0 - 6.5 mg/ml
F/P Molar Ratio	2.5 - 7.5
IF Cells	Rabbit Spleen or Rabbit Blood Lymphocytes
IF-Titer	1:300-1:600
IHC - Tissue	Rabbit Spleen
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:300-1:600
Recommended Retest Period	-----
3 Years	

Specification: PRE.1.ZQ5.10000015780

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.