

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cannabidiol solution - 1.0 mg/mL in methanol, analytical standard, for drug analysis

Product Number:

C6395

MDL:

MFCD00869597

Formula:


C₂₁H₃₀O₂

Formula Weight:

314.46 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Concentration

0.95 - 1.05 mg/ml

Infrared spectrum

Conforms to Structure

Expiration Dated Product

Specification: PRD.2.ZQ5.10000041397

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.