


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
β-Sitosterol - from soybean, ≥96%

Product Number: **S9889**
CAS Number: 83-46-5
MDL: MFCD00003631
Formula: C₂₉H₅₀O
Formula Weight: 414.71 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless
20 mg/mL in CHCL ₃	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 96 %

Specification: PRD.0.ZQ5.10000076227

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.