


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
D-Camphor - ≥97%, FG

Product Number: W223018
CAS Number: 464-49-3
MDL: MFCD00064149
Formula: C₁₀H₁₆O
Formula Weight: 152.23 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Solid
Melting Point	174.0 - 182.0 °C
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.0 %
Specific Rotation	41.00 - 45.00 °
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.5.ZQ5.10000012717

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.