


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,4,7,9-Tetramethyl-5-decyne-4,7-diol, mixture of (±) and meso - 98%

Product Number: 278386
CAS Number: 126-86-3
MDL: MFCD00008942
Formula: C₁₄H₂₆O₂
Formula Weight: 226.36 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White (Powder or Solid) or Colorless to Pale Yellow (Liquid)	
Appearance (Form)	Conforms to Requirements
Powder, Solid or Liquid	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000099905

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.