

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)**Product Specification**

Product Name:

Nitrilotriacetic acid trisodium salt Sigma Grade, ≥98%

Product Number:

**N0253**

CAS Number:

5064-31-3

MDL:


MFCD00064231

Formula:

C<sub>6</sub>H<sub>6</sub>NNa<sub>3</sub>O<sub>6</sub>

Formula Weight:

257.08 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 50 mg/ml, H <sub>2</sub> O	Clear
Proton NMR spectrum	Conforms to Structure
% Water by NMR	≤10.0
Titration calculated as the monohydrate	≥98 %
Sodium (Na) anhydrous	25.0- 27.0 %

**Remarks:**

Specification Date : 01/31/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.