


## Product Specification

**Product Name:**

Kanamycin solution from *Streptomyces kanamyceticus* - 50 - 60 mg/mL in 0.9% NaCl, BioReagent, liquid,  
0.1 µm filtered, suitable for cell culture

**Product Number:****K0254****Storage Temperature:**

2 - 8 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Osmolality	≤ 500 mOs/kg
Concentration Kanamycin	50 - 60 mg/ml
Endotoxin Level at 2ml/L	≤ 0.5 EU/ml
Cell Line	Hep-2 Human Cells
Cell Culture Test of control	≥ 50 %
Expiration Date Period 24 months	-----

Specification: PRE.5.ZQ5.10000014805

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

