

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Acetyl-CoA Carboxylase 1 human – recombinant, expressed in Sf9 cells

Product Number: **A6986**

Storage Temperature: -70 °C

TEST	Specification
------	---------------

mg protein/ml (Bradford)

≥ 0.01

units/μg protein

≥ 20

Enzyme Activity: One unit will cause the conversion of 1 picomole of ATP to ADP per minute at pH 7.4 at 30 deg C

Specification: PRD.6.ZQ5.10000013587

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.