

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Naphthol AS-D Chloroacetate solution

Product Number: 911

MDL: MFCD00042714



TEST

Specification

Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Solution
Appearance (Turbidity)	Clear
Essentially free of particles/precipitate	
Expiration Date	Confirmed
18 months	
Performance	Conforms
When staining with Sigma-Aldrich Procedure 91, sites of activity show	bright red granulation; activity weak or absent in monocytes and lymphocytes

Specification: PRD.2.ZQ5.10000027838

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.