

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Poly(ethylene glycol) linear dendron - NHBoc terminated, PEG Mn 2000, generation 1

Product Number:

901359



Formula:

C₁₁₂H₂₁₈N₂O₅₅

Formula Weight:

2,472.94 g/mol

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

White to Off White

Appearance (Form)

Powder or Crystal or Chunk(s)

1H NMR Spectrum

Conforms to Structure

Molecular Number by HNMR

1000 - 3000

PEG

Recommended Retest Period

2 Years

Specification: PRD.0.ZQ5.10000101983

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.