

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Methyl-1H-pyrazole-4-boronic acid pinacol ester 97%

Product Number:

706078

CAS Number:


936250-20-3

Formula:

C₁₀H₁₇BN₂O₂

Formula Weight:

208.07 g/mol



TEST	Specification
Appearance (Color)	White to Tan
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	55.7- 59.8 %
Nitrogen	13.0- 13.9 %
Purity (GC)	≥96.5 %

Remarks:

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.