

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Bromobutylboronic acid - $\geq 95\%$

Product Number: **691607**
CAS Number: 61632-72-2

Formula: C₄H₁₀BBrO₂
Formula Weight: 180.84 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
Appearance (Form)
Solid, Powder, Crystals, Crystalline Powder, and/or Chunks
Proton NMR Spectrum
Titration with NaOH

White to Tan
Conforms to Requirements

Conforms to Structure
 $\geq 95.0\%$

Specification: PRD.0.ZQ5.10000038695

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.