

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cyclopentyl methyl ether – contains 50 ppm BHT as inhibitor, anhydrous, ≥99.9%

Product Number:

675970

CAS Number:

5614-37-9

MDL:


MFCD08276401

Formula:

C₆H₁₂O

Formula Weight:

100.16 g/mol



TEST

Specification

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Color Test

≤ 10 APHA

Purity (GC)

≥ 99.9 %

BHT as Inhibitor

45 - 90 ppm

Peroxide (ppm)

≤ 50 ppm

Water (by Karl Fischer)

≤ 0.0050 %

Specification: PRD.0.ZQ5.10000020716

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.