

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver acetate – 99.99% trace metals basis

Product Number:

204374

CAS Number:

563-63-3

MDL:


MFCD00012446

Formula:

C₂H₃AgO₂

Formula Weight:

166.91 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Crystalline Powder

Titration with KSCN

63.7 - 65.6 %

% Silver

ICP Major Analysis

Confirmed

Confirms Silver Component

Purity

Meets Requirements

99.99% Based on Trace Metal Analysis

Trace Metal Analysis

≤ 150.0 ppm

Specification: PRD.3.ZQ5.10000020730

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.