


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Methylcatechol - ≥95%

Product Number: M34200
CAS Number: 452-86-8
MDL: MFCD00002205
Formula: C₇H₈O₂
Formula Weight: 124.14 g/mol



TEST

Specification

Appearance (Color)	White to Dark Beige to Dark Brown
Appearance (Form)	Powder or Chunks
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 95.0 %
Solubility (Turbidity) 50mg/ml, Ethanol	Clear to Slightly Hazy
Solubility (Color)	Colorless to Yellow to Orange

Specification: PRD.1.ZQ5.10000010316

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.