


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Pyrrolidinone - 99%

Product Number: P74370
CAS Number: 616-45-5
MDL: MFCD00005270
Formula: C₄H₇NO
Formula Weight: 85.10 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Light Yellow and White to Light Yellow
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %

Specification: PRD.1.ZQ5.10000005046

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.