


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Phenylvinylboronic acid - 95%

Product Number: 571350
CAS Number: 14900-39-1
MDL: MFCD01074578
Formula: C₈H₉BO₂
Formula Weight: 147.97 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
Appearance (Form)
Powder, Crystals, Crystalline Powder, Flakes and/or Chunks
Proton NMR Spectrum
Purity (Titration by NaOH)

Pale Yellow to Dark Yellow
Conforms to Requirements
Conforms to Structure
95.0 - 105.0 %

Specification: PRD.1.ZQ5.10000008931

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.